ETHICON™

SURGIFLO™

HAEMOSTATIC MATRIX

8 mL

en	<i>Instructions for Use</i>	lt	<i>Naudojimo instrukcija</i>
ar	<i>تعليمات الاستخدام</i>	lv	<i>Lietošanas instrukcija</i>
bg	<i>ИНСТРУКЦИИ ЗА УПОТРЕБА</i>	nl	<i>Gebruiksaanwijzing</i>
cs	<i>Návod k použití</i>	no	<i>Bruksanvisning</i>
da	<i>Bruksanvisning</i>	pl	<i>Instrukcja użytkowania</i>
de	<i>Gebrauchsanweisung</i>	pt	<i>Instruções de utilização</i>
el	<i>Οδηγίες χρήσης</i>	ro	<i>Instrucțiuni de utilizare</i>
es	<i>Instrucciones de uso</i>	ru	<i>ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ</i>
et	<i>Kasutusjuhised</i>	sk	<i>Návod na použitie</i>
fi	<i>Käyttöohjeet</i>	sl	<i>Navodila za uporabo</i>
fr	<i>Mode d'emploi</i>	sv	<i>Bruksanvisning</i>
hr	<i>Upute za uporabu</i>	tr	<i>Kullanım Talimatları</i>
hu	<i>Használati utasítás</i>	vi	<i>Hướng dẫn sử dụng</i>
it	<i>Istruzioni per l'uso</i>	zh	<i>使用說明</i>

4548686

REF **MS0010**

English	Page 3
العربية	الصفحة 7
Български	Страница 11
Česky	Strana 16
Dansk	Side 20
Deutsch	Seite 24
Ελληνικά	Σελίδα 28
Español	Página 33
Eesti	Lehekülj 38
Suomi	Sivu 42
Français	page 46
Hrvatski	Stranica 51
Magyar	oldal 55
Italiano	Pagina 59
Lietuvių k.	Puslapis 63
Latviešu	Lappuse 67
Nederlands	Pagina 71
Norsk	Side 75
Polski	Strona 79
Português	Página 84
Română	Pagina 88
Русский	Стр. 92
Slovensky	Strana 97
Slovenščina	Stran 101
Svenska	Sida 105
Türkçe	Sayfa 109
Tiếng Việt	Trang 113
繁體中文	第 117 頁

SURGIFLO™ Haemostatic Matrix

Do not inject into blood vessels.

PRODUCT DESCRIPTION

SURGIFLO™ Haemostatic Matrix (SURGIFLO™) is intended for haemostatic use by applying to a bleeding surface. It is for use only by healthcare professionals who are trained in the surgical procedures and techniques requiring the use of this haemostatic device.

SURGIFLO™ comes in a sterile tray with *all* sterile components for preparing the Flowable Gelatin Matrix:

- A sterile pre-filled blue plunger syringe containing the porcine gelatin matrix that is off-white in appearance
- A sterile empty syringe
- A sterile liquid transfer cup
- A sterile blue flexible applicator tip that is bendable in all directions
- A sterile white applicator tip that can be trimmed to desired length

2 mL of sterile sodium chloride (sterile saline) solution must be added to the Flowable Gelatin Matrix prior to use.

Once the Flowable Gelatin Matrix is mixed with sterile saline, the appropriate applicator tip enclosed in the package may be attached to the syringe for product delivery onto the bleeding site.

The volume following the mixing of the Flowable Gelatin Matrix with the sterile saline is minimum 8 mL.

The clinical benefit to be expected is the control of capillary, venous and arteriolar bleeding when pressure, ligation or other standard surgical methods of control are ineffective or impractical. A Summary of Safety and Clinical Performance for SURGIFLO™ Haemostatic Matrix can be found at the following link in the EUDAMED database:

<https://ec.europa.eu/tools/eudamed>, Basic UDI-DI: 5712123SURGIFLO_w/oThr_HY

ACTIONS

SURGIFLO™ has haemostatic properties. The Flowable Gelatin Matrix provides an environment for platelets to adhere and aggregate within, building on the patient's natural coagulation cascade. When used properly in minimal amounts, SURGIFLO™ is absorbed completely within 4-6 weeks. Absorption is dependent on several factors, including the amount used, degree of saturation with blood or other fluids, and the site of use.

In an animal implantation study, tissue reactions were classified as minimal.

INTENDED USE / INDICATIONS

SURGIFLO™ is indicated for surgical procedures (except ophthalmic) for haemostasis, when control of capillary, venous and arteriolar bleeding by pressure, ligation, and other conventional procedures is ineffective or impractical.

CONTRAINDICATIONS

- Do not inject or compress SURGIFLO™ into blood vessels. Do not use SURGIFLO™ in intravascular compartments because of the risk of thromboembolism, disseminated intravascular coagulation, and increased risk of anaphylactic reaction.
- Do not use SURGIFLO™ in patients with known allergies to porcine gelatin.
- Do not use SURGIFLO™ in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

WARNINGS

- Do not inject or compress SURGIFLO™ into blood vessels as it is for epilesional use only.
- Do not apply SURGIFLO™ in the absence of active blood flow, e.g. while the vessel is clamped or bypassed due to the risk of intravascular clotting from intravascular injection.
- SURGIFLO™ is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for haemostasis. SURGIFLO™ is not intended to be used as a prophylactic haemostatic agent.
- SURGIFLO™ should not be used in the presence of infection. SURGIFLO™ should be used with caution in contaminated areas of the body. If signs of infection or abscess develop where SURGIFLO™ has been positioned, reoperation may be necessary in order to remove or drain the infected material.
- SURGIFLO™ should not be used in instances of pumping arterial haemorrhage. It should not be used where blood or other fluids have pooled or in cases where the point of haemorrhage is submerged. SURGIFLO™ will not act as a tampon or plug in a bleeding site.
- SURGIFLO™ should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. Care should be exercised to avoid overpacking. SURGIFLO™ may swell creating the potential for nerve damage.
- Excess SURGIFLO™ should be removed once haemostasis has been achieved because of the possibility of dislodgement of the device or compression of other nearby anatomic structures.
- The safety and effectiveness of SURGIFLO™ for use in ophthalmic procedures has not been established.
- SURGIFLO™ should not be used for controlling post-partum intrauterine bleeding or menorrhagia.
- The safety and effectiveness of SURGIFLO™ has not been established in children and pregnant or nursing women.

- The blue flexible applicator tip should not be trimmed to avoid exposing internal guidewire.
- The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

PRECAUTIONS

- SURGIFLO™ is for single use only. Do not resterilize. If the product is reused the performance of the product may be deteriorated, and cross contamination may occur which may lead to infection.
- SURGIFLO™ is supplied as a sterile product. Unused open SURGIFLO™ should be discarded. Do not use SURGIFLO™ if sterile barrier package is damaged as sterility may be compromised.
- While packing a cavity for haemostasis is sometimes surgically indicated, SURGIFLO™ should not be used in this manner unless excess product, that is not needed to maintain haemostasis, is removed. When confined in a clot, SURGIFLO™ may swell approximately 20% upon contact with additional fluid.
- Only the minimum amount of SURGIFLO™ needed to achieve haemostasis should be used. Once haemostasis is achieved, any excess SURGIFLO™ should be carefully removed. It is recommended to remove excess of SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.
- SURGIFLO™ should not be used in conjunction with autologous blood salvage circuits. It has been demonstrated that fragments of collagen-based haemostatic agents may pass through 40µm transfusion filters of blood scavenging systems.
- SURGIFLO™ should not be used in conjunction with methylmethacrylate adhesives. Microfibrillar collagen has been reported to reduce the strength of methylmethacrylate adhesives used to attach prosthetic devices to bone surfaces.
- SURGIFLO™ should not be used for the primary treatment of coagulation disorders.
- As with other collagen/gelatin-based haemostatic agents that swell, SURGIFLO™ should be used with caution in neurosurgery. Safe and effective use of SURGIFLO™ in neurosurgery has not been established through randomized, controlled clinical studies.
- Although the safety and effectiveness of the combined use of SURGIFLO™ with other agents has not been evaluated in controlled clinical trials, if in the physician's judgement, concurrent use of other agents is medically advisable, the product literature for that agent should be consulted for complete Prescribing Information.
- The safety and effectiveness of the combined use of SURGIFLO™ with antibiotic solutions or powders has not been established.
- As with other collagen/gelatin-based haemostatic agents, in urological procedures SURGIFLO™ should be removed from the renal pelvis, renal calyces, bladder, urethra or ureters to eliminate the potential foci for calculus formation. The safety and effectiveness of SURGIFLO™ for use in urological procedures has not been established through a randomized clinical study.

SERIOUS INCIDENT REPORTING

Any serious incident that has occurred in relation to the use of SURGIFLO™ or as a result of its use within an EU member state should be reported to Ferrosan Medical Devices A/S through the following email address: complaints@ferrosanmd.com. The user should also report the incident to the national competent authority.

GELATIN-BASED HAEMOSTATIC AGENTS: REPORTED ADVERSE EVENTS

In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based haemostatic agents:

- Gelatin-based haemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.
- Giant cell granulomas have been observed at implant sites when used in the brain.
- Compression of the brain and spinal cord resulting from the accumulation of sterile fluid has been observed.
- Multiple neurologic events were reported when absorbable gelatin-based haemostatic agents were used in laminectomy operations, including cauda equina syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.
- The use of absorbable gelatin-based haemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, cauda equina syndrome, neurogenic bladder, impotence, and paresis.
- The use of absorbable gelatin-based haemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.
- Foreign body reactions, "encapsulation" of fluid, and haematoma have been observed at implant sites.
- Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.
- Toxic shock syndrome was reported in association with the use of absorbable gelatin-based haemostats in nasal surgery.
- Fever, failure of absorption, and hearing loss have been observed when absorbable haemostatic agents were used during tympanoplasty.

ADVERSE REACTIONS REPORTED FROM UNAPPROVED USES

As with other collagen/gelatin-based topical haemostatic agents used for catheter embolization, there is a risk of thromboembolism, pseudoaneurysms and delayed bleeding events, if product is used to obliterate or seal a tract communicating with large vessels.

HOW SUPPLIED

SURGIFLO™ comes in a sterile tray with all sterile components:

- A sterile pre-filled syringe with blue plunger containing the sterile Gelatin Matrix. The syringe is labeled **SURGIFLO™ Hemostatic Matrix**
- A sterile empty syringe
- A sterile liquid transfer cup
- A sterile blue applicator tip that is bendable in all directions
- A sterile white applicator tip that can be trimmed to desired length

The sterile Flowable Gelatin Matrix and all the accessories in the sterile tray are sterilized by gamma irradiation.

SURGIFLO™ also contains Instructions for Use and tracking labels.

STORAGE AND HANDLING

- SURGIFLO™ should be stored dry at controlled temperature (2 °C-25 °C).
- SURGIFLO™ is for single use only.
- The Flowable Gelatin Matrix may be used up to eight (8) hours after mixing with the sterile saline. When the product has been mixed and is ready to use, it is contained in the sterile syringe with the blue plunger that is labeled **SURGIFLO™ Hemostatic Matrix**. The syringe should be kept in the sterile field under room temperature.

DIRECTIONS FOR USE

Before use:


Inspect the sterile barrier package for signs of damage. If the packaging is damaged, opened, or wet, sterility cannot be assured, and the contents should not be used.

Unused opened package of SURGIFLO™ should be discarded since it is not intended for re-use and/or re-sterilization.


Opening the Flowable Gelatin Matrix tray:

Open the outer package and deliver the sterile inner tray to the sterile field using aseptic technique. Once placed in the sterile field, the sterile inner tray may be opened.

Preparing the Flowable Gelatin Matrix with the sterile saline in the sterile field:


1) Draw 2 mL of sterile saline from sterile liquid transfer cup into the empty sterile syringe.

2) Connect syringes

Remove the blue cap from the end of the sterile pre-filled syringe with the blue plunger containing the Flowable Gelatin Matrix. Attach this syringe to the sterile syringe containing the sterile saline.


3) Mix contents of the 2 syringes

Begin mixing by transferring the sterile saline into the sterile pre-filled syringe containing the Flowable Gelatin Matrix. Push the combined material back and forth 6 times until the consistency is even.

Once mixed, the haemostatic matrix should reside completely in the syringe with the blue plunger that is labeled **SURGIFLO™ Hemostatic Matrix**.

Remove the empty syringe and discard.

4) Attach applicator tip

The product is now ready for clinical use.

- The blue flexible applicator tip is bendable in all directions. (Do not trim the flexible applicator tip to avoid exposing the internal guidewire).
- The white applicator tip can be trimmed to desired length. The tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip. The tray can be used to contain the excess piece(s) for discarding.

Do not inject SURGIFLO™ into blood vessels. See the Contraindications, Warnings and Precautions.

For open procedures:

- a. Identify the source of bleeding.
- b. Deliver SURGIFLO™ to the source of bleeding. SURGIFLO™ can be used with or without one of the applicator tips attached to the syringe labeled **SURGIFLO™ Hemostatic Matrix**. Apply sufficient SURGIFLO™ to cover the entire bleeding surface.

- c. For tissue defects (cavities, divots, or craters), apply SURGIFLO™ at the deepest part of the lesion, and continue applying material as the syringe (or applicator tip) is withdrawn from the lesion.
- d. Apply a sterile saline moistened gauze over the SURGIFLO™ to ensure the material remains in contact with the bleeding tissue.
- e. After 1-2 minutes, lift and remove the gauze and inspect the bleeding site.
- f. When bleeding has ceased, remove the gauze and excess SURGIFLO™ should be removed. It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.
- g. In cases of persistent bleeding indicated by saturation and bleeding through the material, it is important that the surgeon performs a re-assessment of the bleeding site to determine the appropriate management to achieve haemostasis.

If it is evaluated to be clinically appropriate to re-apply SURGIFLO™, the product may be re-applied to the bleeding site following the steps in point d., e., and f.

It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.

For endoscopic sinus surgery and epistaxis:

- a. Deliver SURGIFLO™ to the source of bleeding using the selected applicator tip attached to the **SURGIFLO™ Hemostatic Matrix** syringe.
- b. Apply sufficient SURGIFLO™ to cover the entire bleeding surface.
- c. Using forceps or an appropriate instrument, carefully layer a sterile saline moistened gauze over the SURGIFLO™ for 1-2 minutes to ensure the material remains in contact with the bleeding tissue.
- d. After 1-2 minutes, lift and remove the gauze and inspect the bleeding site.
- e. When bleeding has ceased, remove the gauze and excess SURGIFLO™ should be removed. It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.
- f. In cases of persistent bleeding indicated by saturation and bleeding through the material, it is important that the surgeon performs a re-assessment of the bleeding site to determine the appropriate management to achieve haemostasis.

If it is evaluated to be clinically appropriate to re-apply SURGIFLO™, the product may be re-applied to the bleeding site following the steps in point c., d., and e.


It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.

- g. Use of nasal packing is not necessary when satisfactory haemostasis is achieved.

DISPOSAL OF SURGIFLO™ AFTER USE

Dispose any remaining SURGIFLO™ Flowable Gelatin Matrix, the accessory components, and the packaging according to your facility's policies and procedures regarding biohazardous materials and waste.

SYMBOLS USED ON LABELING

	Consult instructions for use		Medical Device
	Not made with natural rubber latex		Contains biological material of animal origin
	Contains porcine gelatin as the main ingredient		Distributor
	Do not inject into blood vessels		Manufacturer
	Single sterile barrier system with protective packaging inside; device inside this packaging is sterile		CE mark and identification number of notified body
	Sterile; the medical device is supplied sterile		
	Sterilized using irradiation		Reorder number
	Do not use if package is damaged and consult Instructions for Use		Batch number
	Do not re-use		Date of manufacture year, month and day
	Do not re-sterilize		Use by year, month and day
	Temperature limit		Designates the packaging material to which it is applied is recyclable. Recycling programs may not exist in your area

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© Ethicon, Inc. 2022
Ethicon, Inc.