Johnson&Johnson

# SURGIFLO™

Haemostatic Matrix Kit

with Thrombin

**8 mL**

**en** *Instructions for Use*

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# SURGIFLO™ Haemostatic Matrix Kit with Thrombin

**Do not inject into blood vessels.**

## PRODUCT DESCRIPTION

SURGIFLO™ Haemostatic Matrix Kit with Thrombin (SURGIFLO™) is intended for haemostatic use by applying to a bleeding surface. It is for use only by healthcare professionals who are trained in the surgical procedures and techniques requiring the use of this haemostatic device.

The Kit contains:

1. A sterile tray (SURGIFLO™ Haemostatic Matrix) with *all* sterile components to prepare the Gelatin Matrix
2. A sterile tray (Thrombin Constituents) with *all* surface sterilized constituents to prepare the Thrombin Solution

SURGIFLO™ Haemostatic Matrix comes in a tray with *all* sterile components:

- A sterile syringe, labeled SURGIFLO™, containing the porcine Gelatin Matrix that is off-white in appearance,
- A sterile blue flexible applicator tip that is bendable in all directions, and
- A sterile white straight applicator tip that can be trimmed to desired length.

The surface sterilized constituents to prepare the Thrombin Solution include:

- A Thrombin vial containing 2000 International Units (IU) of sterile lyophilized human thrombin,
- A needle-free syringe containing 2 mL of sterile water (sWFI) for Thrombin reconstitution. This syringe has a blue flange that is rotatable allowing flexibility to position the syringe for individual preference, if needed,
- A vial adapter.

Thrombin should be reconstituted using the vial adapter and the needle-free syringe with sterile water (sWFI) for Thrombin reconstitution.

The Thrombin Solution must be added to the Gelatin Matrix prior to use. The volume following the mixing of the Gelatin Matrix with the Thrombin Solution is minimum 8 mL. Once the Gelatin Matrix is mixed with the Thrombin Solution, the appropriate applicator tip enclosed in the package may be attached to the syringe for product delivery onto the bleeding site.

The clinical benefit to be expected is the control of bleeding ranging from oozing to spurting when pressure, ligature or other standard surgical methods of control are ineffective or impractical. A Summary of Safety and Clinical Performance for SURGIFLO™ Haemostatic Matrix Kit with Thrombin can be found at the following link in the EUDAMED database: <https://ec.europa.eu/tools/eudamed>, Basic UDI-DI: 5712123SURGIFLO\_w/Thr\_G9

## ACTIONS

SURGIFLO™ has haemostatic properties. The Gelatin Matrix provides an environment for platelets to adhere and aggregate within, building on the patient's natural coagulation cascade.

The patient's endogenous thrombin is activated, and the patient's thrombin converts fibrinogen into an insoluble fibrin clot. The thrombin constituent of SURGIFLO™ provides an ancillary effect to the innate haemostatic property of the Gelatin Matrix.

When used properly in minimal amounts, SURGIFLO™ is absorbed completely within 4-6 weeks. Absorption is dependent on several factors, including the amount used, degree of saturation with blood or other fluids, and the site of use.

In an animal implantation study, tissue reactions were classified as minimal.

## INTENDED USE / INDICATIONS

SURGIFLO™ is indicated in surgical procedures (other than ophthalmic) as an adjunct to haemostasis when control of bleeding ranging from oozing to spurting by ligature or other conventional methods is ineffective or impractical.

## CONTRAINDICATIONS

- Do not inject or compress SURGIFLO™ into blood vessels. Do not use SURGIFLO™ in intravascular compartments because of the risk of thromboembolism, disseminated intravascular coagulation, and increased risk of anaphylactic reaction.
- Do not use SURGIFLO™ in patients known to have anaphylactic or severe systemic reaction to human blood products.
- Do not use SURGIFLO™ in patients with known allergies to porcine gelatin.
- Do not use SURGIFLO™ in closure of skin incisions because it may interfere with the healing of skin edges. This interference is due to mechanical interposition of gelatin and is not secondary to intrinsic interference with wound healing.

## WARNINGS

- Do not inject or compress SURGIFLO™ into blood vessels as it is for epilesional use only.
- Do not apply SURGIFLO™ in the absence of active blood flow, e.g. while the vessel is clamped or bypassed due to the risk of intravascular clotting from intravascular injection.
- SURGIFLO™ contains thrombin made from human plasma. Products made from human plasma may carry a risk of transmitting infectious agents, such as viruses, and theoretically, the Creutzfeldt-Jakob disease (CJD) agent. The risk of transmitting an infectious agent has been reduced by screening plasma donors for prior exposure to certain viruses, by testing for the presence of certain current virus infections, and by inactivating and removing certain viruses. Despite these measures, such products can still potentially transmit disease. There is also the possibility that unknown infectious agents may be present in such products. The physician should discuss the risks and benefits of this product with the patient. The measures taken are considered effective for enveloped viruses such as HIV, HCV, HBV, and for the non-enveloped virus HAV. The measures taken may be of limited value against non-enveloped viruses such as parvovirus B19. Parvovirus B19 infection may be serious for pregnant women (foetal infection) and for individuals with immunodeficiency or increased erythropoiesis (e.g. haemolytic anaemia).
- SURGIFLO™ is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for haemostasis. SURGIFLO™ is not intended to be used as a prophylactic haemostatic agent.
- SURGIFLO™ should not be used in the presence of infection. SURGIFLO™ should be used with caution in contaminated areas of the body. If signs of infection or abscess develop where SURGIFLO™ has been positioned, reoperation may be necessary in order to remove or drain the infected material.
- SURGIFLO™ should not be used in instances of pumping arterial haemorrhage. It should not be used where blood or other fluids have pooled or in cases where the point of haemorrhage is submerged. SURGIFLO™ will not act as a tampon or plug in a bleeding site.
- SURGIFLO™ should be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. Care should be exercised to avoid overpacking. SURGIFLO™ may swell creating the potential for nerve damage.

- Excess SURGIFLO™ should be removed once haemostasis has been achieved because of the possibility of dislodgement of the device or compression of other nearby anatomic structures.
- The safety and effectiveness of SURGIFLO™ for use in ophthalmic procedures has not been established.
- SURGIFLO™ should not be used for controlling post-partum intrauterine bleeding or menorrhagia.
- The safety and effectiveness of SURGIFLO™ has not been established in children and pregnant or nursing women.
- The blue flexible applicator tip should not be trimmed to avoid exposing the internal guidewire.
- The white straight applicator tip should be trimmed away from the surgical area. Cut a square angle to avoid creating a sharp tip.

## PRECAUTIONS

- SURGIFLO™ is for single use only. Do not resterilize. If the product is reused the performance of the product may be deteriorated, and cross contamination may occur which may lead to infection.
- SURGIFLO™ is supplied as a sterile product. Unused open SURGIFLO™ should be discarded. Do not use SURGIFLO™ if sterile barrier package is damaged as sterility may be compromised.
- While packing a cavity for haemostasis is sometimes surgically indicated, SURGIFLO™ should not be used in this manner unless excess product, that is not needed to maintain haemostasis, is removed. When confined in a clot, SURGIFLO™ may swell approximately 20% upon contact with additional fluid.
- Only the minimum amount of SURGIFLO™ needed to achieve haemostasis should be used. Once haemostasis is achieved, any excess SURGIFLO™ should be carefully removed. It is recommended to remove excess of SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.
- SURGIFLO™ should not be used in conjunction with autologous blood salvage circuits. It has been demonstrated that fragments of collagen-based haemostatic agents may pass through 40µm transfusion filters of blood scavenging systems.
- SURGIFLO™ should not be used in conjunction with methylmethacrylate adhesives. Microfibrillar collagen has been reported to reduce the strength of methylmethacrylate adhesives used to attach prosthetic devices to bone surfaces.
- Similar to comparable products containing thrombin, the Thrombin Solution may be denatured after exposure to solutions containing alcohol, iodine or heavy metals (e.g. antiseptic solutions). Such substances should be removed to the greatest possible extent before applying the product.
- SURGIFLO™ should not be used for the primary treatment of coagulation disorders.
- As with other collagen/gelatin based-haemostatic agents, in urological procedures SURGIFLO™ should not be left in the renal pelvis, renal calyces, bladder, urethra or ureters to eliminate the potential foci for calculus formation. The safety and effectiveness of SURGIFLO™ for use in urological procedures has not been established through a randomized clinical study.
- As with other collagen/gelatin based-haemostatic agents that swell, SURGIFLO™ should be used with caution in neurosurgery. Safe and effective use of SURGIFLO™ in neurosurgery has not been established through randomized, controlled clinical studies.
- Although the safety and effectiveness of the combined use of SURGIFLO™ with other agents have not been evaluated in controlled clinical trials, if in the physician's judgement, concurrent use of other agents is medically advisable, the product literature for that agent should be consulted for complete Prescribing Information.
- The safety and effectiveness of the combined use of SURGIFLO™ with antibiotic solutions or powders has not been established.

## SERIOUS INCIDENT REPORTING

Any serious incident that has occurred in relation to the use of SURGIFLO™ or as a result of its use within an EU member state should be reported to Ferrosan Medical Devices A/S through the following email address: [complaints@ferrosanmd.com](mailto:complaints@ferrosanmd.com). The user should also report the incident to the national competent authority.

## GELATIN-BASED HAEMOSTATIC AGENTS: REPORTED ADVERSE EVENTS

In general, the following adverse events have been reported with the use of absorbable porcine gelatin-based haemostatic agents:

- Gelatin-based haemostatic agents may serve as a nidus for infection and abscess formation and have been reported to potentiate bacterial growth.
- Giant cell granulomas have been observed at implant sites when used in the brain.
- Compression of the brain and spinal cord resulting from the accumulation of sterile fluid has been observed.
- Multiple neurologic events were reported when absorbable gelatin-based haemostatic agents were used in laminectomy operations, including, cauda equina syndrome, spinal stenosis, meningitis, arachnoiditis, headaches, paresthesias, pain, bladder and bowel dysfunction, and impotence.
- The use of absorbable gelatin-based haemostatic agents during the repair of dural defects associated with laminectomy and craniotomy operations has been associated with fever, infection, leg paresthesias, neck and back pain, bladder and bowel incontinence, cauda equina syndrome, neurogenic bladder, impotence, and paresis.
- The use of absorbable gelatin-based haemostatic agents has been associated with paralysis, due to device migration into foramina in the bone around the spinal cord, and blindness due to device migration in the orbit of the eye, during lobectomy, laminectomy, and repair of a frontal skull fracture and lacerated lobe.
- Foreign body reactions, "encapsulation" of fluid, and haematoma have been observed at implant sites.
- Excessive fibrosis and prolonged fixation of a tendon have been reported when absorbable gelatin-based sponges were used in severed tendon repair.
- Toxic shock syndrome was reported in association with the use of absorbable gelatin-based haemostats in nasal surgery.
- Fever, failure of absorption, and hearing loss have been observed when absorbable haemostatic agents were used during tympanoplasty.

## ADVERSE REACTIONS TO HUMAN THROMBIN

As with any other plasma derivatives, hypersensitivity or allergic reactions may occur in rare cases. In isolated cases, these reactions may progress to severe anaphylaxis. Other adverse events reported in a clinical trial were abnormal laboratory tests (prolonged activated partial thromboplastin time, prolonged prothrombin time, increased INR, decreased lymphocyte count, increased neutrophil count) and haematoma.

## ADVERSE REACTIONS TO GELATIN-BASED HAEMOSTATIC AGENTS WITH THROMBIN

Adhesion formation and small bowel obstruction are well-known and common complications following abdominal and gynaecological surgeries.

Adverse events such as these, as well as inflammation and foreign body reaction including giant cell granulomas, have been reported within abdominal and gynaecological surgeries in which gelatin-based and collagen-based haemostatic agents with thrombin have been used in excessive amount.

Therefore, as with other gelatin haemostatic agents, only the minimum amount of SURGIFLO™ needed to achieve haemostasis should be used. Once haemostasis is achieved, any excess SURGIFLO™ should be carefully removed.

## ADVERSE REACTIONS REPORTED FROM UNAPPROVED USES

As with other collagen/gelatin-based topical haemostatic agents used for catheter embolization, there is a risk of thromboembolism, pseudoaneurysms and delayed bleeding events, if product is used to obliterate or seal a tract communicating with large vessels.

## HOW SUPPLIED

SURGIFLO™ consists of:

1. A sterile tray (SURGIFLO™ Haemostatic Matrix) with *all* sterile components to prepare the Gelatin Matrix as described in the table below
2. A sterile tray (Thrombin Constituents) with *all* surface sterilized constituents to prepare the Thrombin Solution as described in the table below

SURGFLO™ is provided in the configuration shown in the table below.

SURGFLO™ Haemostatic Matrix Kit with Thrombin	
SURGFLO™ Haemostatic Matrix	Thrombin Constituents
<ul style="list-style-type: none"> <li>• A sterile pre-filled syringe, labeled SURGFLO™, containing the porcine Gelatin Matrix</li> <li>• A sterile blue flexible applicator tip</li> <li>• A sterile white straight applicator tip</li> </ul>	<ul style="list-style-type: none"> <li>• A Thrombin vial containing 2000 International Units (IU) of sterile lyophilized human thrombin for reconstitution</li> <li>• A needle-free syringe containing 2 mL of sterile water (sWFI) for Thrombin reconstitution. This syringe has a rotatable blue flange and is labeled SURGFLO™ Haemostatic Matrix</li> <li>• A vial adapter</li> </ul>

**The sterile Gelatin Matrix and the accessories:**

The tray is sterilized by gamma irradiation.

**The Thrombin kit constituents:**

The tray is surface sterilized with ethylene oxide.

- The lyophilized thrombin (human) is sterilized through an autoclaved sterile filter cartridge.
- The sterile water (sWFI) for Thrombin reconstitution in the needle-free syringe is steam sterilized.

SURGFLO™ contains Instructions for Use. Tracking labels are provided in the kit to record name and batch number of the product to link use to the patient record. It is strongly recommended that every time that SURGFLO™ is administered to a patient, the name and batch number of the product are recorded in order to maintain a link between the patient and the batch of the product.

**STORAGE AND HANDLING**

- SURGFLO™ should be stored dry at controlled temperature (2 °C-25 °C).
- SURGFLO™ is for single use only.
- The Thrombin vial should be kept away from light.
- The Thrombin Solution should be used together with the Gelatin Matrix and used only as indicated.
- SURGFLO™ may be used up to eight (8) hours after mixing with the Thrombin Solution. When the product has been mixed and is ready to use, it is contained in the sterile syringe with the blue flange that is labeled **SURGFLO™ Haemostatic Matrix**. The syringe should be kept in the sterile field under room temperature.

**DIRECTIONS FOR USE**


**Before use:**

Inspect the sterile barrier package for signs of damage. If the packaging is damaged, opened, or wet, sterility cannot be assured, and the contents should not be used. Unused opened packages of SURGFLO™ should be discarded since they are not intended for re-use and/or resterilization.


**Opening the tray with Gelatin Matrix and the tray with Thrombin kit constituents:**

Open the outer packages and deliver the sterile inner trays to the sterile field using aseptic technique. Once placed in the sterile field, the sterile inner tray may be opened.


**Preparing the Thrombin Solution inside the sterile field:**


1. Flip off the cap from the Thrombin vial, leaving the aluminum ring and the rubber stopper in place.


2. Hold or place the Thrombin vial on a flat surface and use the vial adapter to pierce the rubber stopper. Ensure that the vial adapter is fully inserted into the Thrombin vial.

3. Remove the dark blue cap from the end of the needle-free syringe with blue flange containing sterile water (sWFI) for Thrombin reconstitution. Connect and screw the syringe on the vial adapter. Transfer the entire volume of sterile water (sWFI) for Thrombin reconstitution into the Thrombin vial.


4. Without removing the syringe, gently swirl the Thrombin vial until the Thrombin is completely dissolved. This may take about four (4) swirls.


5. Withdraw the sterile Thrombin Solution into the needle-free syringe with the blue flange.

**Preparing the Gelatin Matrix with the Thrombin Solution inside the sterile field:**

**1. Connect syringes**


Remove the white cap from the end of the sterile pre-filled syringe containing the Gelatin Matrix. Attach this syringe to the sterile syringe with the blue flange containing the sterile Thrombin Solution.

**2. Mix contents of the 2 syringes**

Transfer the sterile Thrombin Solution into the pre-filled syringe containing the Gelatin Matrix. Begin mixing by pushing the combined material back and forth for a total of 6 times. Transfer the mixed haemostatic material to the syringe with the blue flange labelled **SURGIFLO™ Haemostatic Matrix**.

Remove the empty syringe and discard.

**3. Attach applicator tip**

The sterile blue flexible applicator or the sterile white straight applicator tip that can be trimmed may be attached to the syringe labeled **SURGIFLO™ Haemostatic Matrix**.

The product is now ready for clinical use.

Do not inject SURGIFLO™ into blood vessels. See the Contraindications, Warnings and Precautions.

**Application of SURGIFLO™:**

**For open procedures:**

- Identify the source of bleeding.
- Deliver SURGIFLO™ to the source of bleeding. SURGIFLO™ can be used with or without one of the applicator tips attached to the syringe labeled **SURGIFLO™ Haemostatic Matrix**. Apply sufficient SURGIFLO™ to cover the entire bleeding surface.
- For tissue defects (cavities, divots or craters) apply SURGIFLO™ at the deepest part of the lesion and continue applying material as the syringe (or applicator tip) is withdrawn from the lesion.
- Apply a sterile saline moistened gauze over the SURGIFLO™ to ensure the material remains in contact with the bleeding tissue.
- After 1-2 minutes, lift and remove the gauze and inspect the bleeding site.
- When bleeding has ceased, remove the gauze and excess SURGIFLO™ should be removed. It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.
- In cases of persistent bleeding indicated by saturation and bleeding through the material, it is important that the surgeon performs a re-assessment of the bleeding site to determine the appropriate management to achieve haemostasis.

If it is evaluated to be clinically appropriate to re-apply SURGIFLO™, the product may be re-applied to the bleeding site following the steps in point d., e., and f.

It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.

**For endoscopic sinus surgery and epistaxis:**

- Deliver SURGIFLO™ to the source of bleeding using the selected applicator tip attached to the **SURGIFLO™ Haemostatic Matrix** syringe.
- Apply sufficient SURGIFLO™ to cover the entire bleeding surface.
- Using forceps or an appropriate instrument, carefully layer a sterile saline moistened gauze over the SURGIFLO™ to ensure the material remains in contact with the bleeding tissue.
- After 1-2 minutes, lift and remove the gauze and inspect the bleeding site.
- When bleeding has ceased, remove the gauze and excess SURGIFLO™ should be removed. It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.
- In cases of persistent bleeding indicated by saturation and bleeding through the material, it is important that the surgeon performs a re-assessment of the bleeding site to determine the appropriate management to achieve haemostasis.

If it is evaluated to be clinically appropriate to re-apply SURGIFLO™, the product may be re-applied to the bleeding site following the steps in point c. d., and e.


It is recommended to remove excess SURGIFLO™ with irrigation and aspiration once haemostasis is achieved, without disturbing the clot.

- Use of nasal packing is not necessary when satisfactory haemostasis is achieved.

## DISPOSAL OF SURGIFLO™ AFTER USE

Dispose any remaining SURGIFLO™ Gelatin Matrix, the accessory components, the thrombin constituent parts (the thrombin vial is made of glass) and the packaging according to your facility's policies and procedures regarding biohazardous materials and waste.

## SYMBOLS USED ON LABELING

 Consult instructions for use	 Keep away from sunlight
 Not made with natural rubber latex	 Medical Device
 Contains porcine gelatin as the main ingredient	 Contains biological material of animal origin
 Do not inject into blood vessels	 Contains human blood or plasma derivatives
 Single sterile barrier system with protective packaging inside; device inside this packaging is sterile	 Contains a medicinal substance
 Sterile; the medical device is supplied sterile	 Distributor
 Sterilized using irradiation	 Manufacturer
 Sterilized using ethylene oxide	 CE mark and identification number of notified body 2460
 Sterilized using aseptic processing techniques	 Reorder number
 Sterilized using steam or dry heat	 Batch number
 Do not use if package is damaged and consult Instructions for Use	 Date of manufacture year, month and day
 Do not re-use	 Use by year, month and day
 Do not re-sterilize	 Indicates the number of pieces in the package
 Temperature limit 2°C - 25°C	 Designates the packaging material to which it is applied is recyclable. Recycling programs may not exist in your area

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© Johnson & Johnson and its affiliates 2025
Ethicon, Inc.